# Pulse Oximeter User Guide ACCURO



(E

**Charmcare Co., Ltd.** 

# **Table Of Contents**

| Part 1 | Before You Begin                              | 5  |
|--------|---|----|
|        | Overview                                      | 5  |
|        | Intended Audience                             | 5  |
|        | Safety Information, Warning, Caution and Note | 5  |
|        | Safety Symbol Definitions                     | 5  |
|        | Warning                                       | 6  |
|        | Caution                                       | 7  |
|        | Proper Environments for the Product           | 9  |
|        | Warranty Period                               | 10 |
|        | Battery Replacement                           | 10 |
|        | Contact Us                                    | 10 |
|        | Electrical Safety                             | 11 |
|        | Maintenance and Cleaning                      | 12 |
|        | Classifications                               | 13 |
|        | Description of Product & Label Symbols        | 14 |
| Part 2 | Product Summary                               | 15 |
|        | Intended Use                                  | 15 |
|        | List of Parts                                 | 18 |
|        | Nomenclature of the Parts                     | 19 |
|        | Unit  | 19 |
| Part 3 | Product Installation                          | 22 |
|        | Caution When Installing                       | 22 |
|        | Connecting the AC Power                       | 22 |
|        | Installing the IV Pole                        | 24 |
| Part 4 | Using The Product                             | 25 |
|        | Preparations before Use                       | 25 |
|        | The Screen                                    | 26 |
|        | SpO2/PR Data Area                             | 27 |

| Pulse Oximeter ACCURO |                                 | www.charmcare.com |
|-----------------------|---------------------------------|-------------------|
|                       | System Status                   | 28                |
|                       | Post-Use Storage and Management | 28                |
|                       | Power Off                       | 29                |
|                       |                                 |                   |
| Part 5                | Measuring the SpO2              | 30                |
|                       | Attaching the SpO2 Probe        | 30                |
|                       | SpO2 Sensor Port                | 31                |
|                       | Measuring the SpO2              | 32                |
|                       |                                 |                   |
| Part 6                | Switch Screen Mode              | 33                |
|                       | Switch Screen Mode              | 33                |
|                       |                                 |                   |
| Part 7                | How to Use Menu                 | 34                |
|                       | Entering Main Menu              | 34                |
|                       | Setup Configuration             | 34                |
|                       | Setup Current Time              | 35                |
|                       | Setup Alarm                     | 36                |
|                       | Setup Patient                   | 37                |
|                       | Entering Maintenance Menu       | 37                |
|                       | Setup Maintenance               | 38                |
|                       |                                 |                   |
| Part 8                | Save And View TREND Data        | 39                |
|                       | Saving Measured Data            | 39                |
|                       | Setup Trend                     | 39                |
|                       | Short Trend View                | 40                |
|                       | Tabular Trend View              | 40                |
|                       | Graphic Trend View              | 41                |
|                       |                                 |                   |
| Part 9                | Alarm                           | 42                |
|                       | Alarm Categories                | 42                |
|                       | Alarm Levels                    | 42                |
|                       | Alarm Indicators                | 42                |
|                       | Priority Signal Timing          | 44                |
|                       | Alarm Conditions                | 45                |

Rev. 2 -3- OP-EN-02

| Pulse Oximeter ACCURO www.charmcare.com |   |    |  |
|---|---|----|--|
|   | Alarm Conditions Delay                        | 45 |  |
|   | •   |    |  |
|   | Alarm Status                                  | 46 |  |
|   | Verifying Visual and Audible Alarm Indication | 48 |  |
|   |   |    |  |
| Part 10                                 | Basic Troubleshooting                         | 49 |  |
|   |   |    |  |
| Part 11                                 | Product Specification                         | 50 |  |
|   | Product Specification                         | 50 |  |
|   | Manufacturer's Declaration                    | 52 |  |
|   | Clinical Studies                              | 55 |  |
|   |   |    |  |
| Attachme                                | ent 1 Factory Set                             | 57 |  |
|   | Factory Set                                   | 57 |  |
|   | Maintenance Set                               | 58 |  |
|   | Factory Default of Patient Type               | 59 |  |
|   |   |    |  |
| Attachme                                | ent 2 Alarm Message                           | 60 |  |
|   |   |    |  |

Rev. 2 -4- OP-EN-02

## Part 1

# Before You Begin

#### 1.1 Overview

This manual contains information for collecting patient oxygen saturation data while operating the ACCURO.

#### 1.2 Intended Audience

This manual provides information to health-care professionals acting as caregivers for operation and user of the monitoring system. Before use, carefully read this manual, accessory Directions for Use, and all precautionary information and specifications.

### 1.3 Safety Information, Warning, Caution and Note

This is contains safety information requiring users to exercise appropriate caution while using the monitoring system.

#### 1.3.1 Safety Symbol Definitions



#### Warning

"Warning" is used to refer to factors, which, when ignored, may result in severe and/or fatal injuries and property damage.



#### Caution

"Caution" is used to refer to factors, which, when ignored, may result in moderate, but non-life-threatening, injuries.



#### Note

"Note" is used to highlight factors that are not dangerous, but should be paid close attention to during installation, use, and maintenance.

Rev. 2 -5- OP-EN-02

#### 

- Do not use the battery with other manufacturer's batteries.
- Do not use any monitoring system or pulse oximeter cables, sensors, or connectors that appear damaged.
- Do not use any monitoring system, sensor, cable, or connector that appears damaged. Remove any damaged equipment from service for inspection by a qualified service technician.
- No modification of this equipment is allowed.
   Do not lift by the sensor or interface cable. The cable may disconnect, potentially dropping the monitoring system on a patient or damaging surface.
- When installing the AC power cord, ensure the cord is carefully positioned to prevent tripping and entanglement.
- Do not spray, pour, or spill any liquid on the monitoring system, its accessories, connectors, switches, or openings in the chassis, since this may cause damage to the monitoring system.
- To ensure accurate performance and prevent device failure, do not subject to extreme moisture, such as direct exposure to rain. Such exposure may cause inaccurate performance or device failure.
- Use only when connected to a grounded outlet to avoid electric shock.
- Do not pause or disable audible alarms or decrease the volume of the audible alarm if patient safety could be compromised. Do not dim or disable visual alarms if patient safety could be compromised.
- Use of portable phones or other radio frequency (RF) emitting equipment near the system may cause unexpected or adverse operation.
- Operator shall not touch the battery compartment and the patient simultaneously.
- The measured values of the monitoring system can be affected by patient conditions, excessive patient movement, sensors, environmental conditions, and nearby electromagnetic external conditions.
- Ensure the monitoring system is clear of any obstructions that prevent awareness of visual or audible alarms. Failure to do so may result in inadvertently missing a visual alarm or an inaudible alarm tone.
- The use of accessories, transducers and cable other than those specified may result in increased emissions or decreased immunity performance of the equipment.
- The equipment or system should not be used adjacent to, or stacked with, other equipment. If adjacent or stacked use is necessary, the equipment or system should be tested to verify normal operation in the configuration in which it is being used.

Rev. 2 -6- OP-EN-02

## 1.3.3 **(** Caution

- Accessory equipment connected to the monitoring system's data interface must be certified
  according to IEC Standard 60950-1 for data-processing equipment. All combinations of
  equipment must be in compliance with IEC Standard 60601-1:2012 Requirements for Medical
  Electrical Systems. Anyone who connects additional equipment to the signal input or signal
  output port configures a medical system and is therefore responsible for ensuring the system
  complies with the requirements of IEC Standard 60601-1:2012 and IEC Standard 60601-12:2007.
- When connecting the monitoring system to any instrument, verify proper operation before clinical use.
- Both the monitoring system and the instrument connected to it must utilize a grounded outlet.
- Any equipment connected to the data interface must be certified according to the latest IEC/EN 60950-1 standard for data-processing equipment, the latest IEC/EN 60601-1 standard for electromedical equipment, or the latest IEC/EN safety standards relevant to that equipment.
- All combinations of equipment must be in compliance with Requirements for Medical Electrical Systems IEC Standard 60601-1:2012.
- Anyone who connects equipment to the data interface is configuring a medical system and, therefore, is responsible for ensuring the system complies with the Requirements for Medical Electrical Systems IEC/EN Standard 60601-1 and the electromagnetic compatibility IEC/EN Standard 60601-1-2:2007.
- Accuracy may degrade if it is connected to secondary I/O devices when the equipment is not connected to earth reference.
- Use only approved sensors and interface cables when connecting to the sensor port. Connecting
  any other cable or sensor influences the accuracy of sensor data, which may lead to adverse
  results.
- The top of the monitoring system screen indicates the sensor type when connecting a recommended sensor to the monitoring system or when the monitoring system completes POST with an attached sensor.
- If the pulse beep tone does not sound with each pulse, the pulse beep volume is set to zero, the speaker is malfunctioning, or the signal is corrupt. Reset the device.
- Perform the following checks every 12 months.
- Inspect the monitoring system for mechanical and functional damage or deterioration.
- Inspect the internal fuses for proper value and rating.
- Ensure all user interface items, cables, and accessories function normally.
- Have a qualified service technician replace the battery at least every 6 months.
- The battery is recyclable. Do not dispose of the battery by placing it in the regular trash.

Rev. 2 -7- OP-EN-02

- Dispose of the battery in accordance with local guidelines and regulations or contact Charmcare to arrange for disposal.
- Inspect the safety relevant labels for legibility. Contact Charmcare or a local representative, if labels are damaged or illegible.
- Instructions provided on preventive inspection, calibration, maintenance and its frequency
- Information provided for safe performance of routine maintenance necessary to ensure continued safe use of ME EQUIPMENT
- Instructions provided to ensure adequate maintenance of ME EQUIPMENT containing rechargeable batteries to be maintained by anyone other than SERVICE PERSONNEL
- Magnetic and electrical fields are capable of interfering with the proper performance of the device.
- X-ray equipment or MRI devices are possible sources of interference as they may emit higher levels of electromagnetic radiation.

Rev. 2 -8- OP-EN-02

#### 1.3.4 Proper Environments For The Product

Do not use or store the product in the following conditions:



#### **Exposure to Humidity and/or Moisture**

Do not use the product with wet hands.



Do not store or place the product in direct sunlight.





Do not store or place the product in areas where temperature changes are drastic.

Use the product between temperatures of 5°C and 35°C and at a humidity between 30% and 85%.



Do not place the product near heat sources.



Do not store or place the product in very humid areas or areas where air circulation is a problem.



Do not subject product to severe shock and or vibrations.





Do not store or place the product in places where product is exposed to chemicals or flammable gas.



Keep product free of dust and debris, particularly metallic objects.





Do not attempt to disassemble the product yourself. Charmcare is not liable for any problems that may occur should you attempt to



Do not connect the power during installation. This may damage the product.





Grip the plug when unplugging the product from a power source.

Rev. 2 -9- OP-EN-02

#### 1.3.5 Warranty Period

- To obtain information about a warranty, if any, for this product, contact your local Charm-care representatives.
- This product has been manufactured and inspected following the strict quality assurance guidelines of Charmcare.
- Refer to the Economic Planning Board's "Regulations Regarding Consumer Compensation" for more information on conditions for product repairs and exchanges.
- Product malfunctions occurring from regular use shall be repaired for free at the Charmcare service center during the term of the warranty period.
- During the term of the warranty period, report all problems with the product to Charmcare by including the model no., the device no., date of purchase and a detailed description of the error.
- Manufacturer and/or the store where the product was purchased do not assume any responsibility for any and/or all problems resulting from improper use or improper storage of the product.

#### 1.3.6 Battery Replacement

- If the monitor has not been used for a long period of time, the batter will need charging.
   To charge the battery, contact your local Charmcare representatives.
- Charmcare recommends that the Li-ion rechargeable battery be replaced at 6 months intervals.
   Refer to the monitor service manual for batter replacement and general service instruction.
   Follow local governing ordinances and recycling plans regarding disposal or recycling of the battery and other device components.
- Remove the battery if the equipment is not likely to be used for some time.
- When incorrect battery replacement by inadequately trained personnel would result in an unacceptable risk(e.g. excessive temperatures, fire or explosion)

#### 1.3.7 Contact Us

- Please contact us for better service and products.
- Charmcare service is always open.

| Need to purchase products | Charmcare Co., Ltd.                                      |  |
|---------------------------|--|--|
| or parts?                 | (Gasan-dong, Woolim Lions2-cha)714, 2, Gasandigital1-ro, |  |
| Need service or repairs?  | Geumcheon-gu, Seoul, Korea                               |  |
| Need technical advice?    | Tel: +82-2-862-5052 Fax: +82-2-862-5065                  |  |
| Website                   | http://www.charmcare.com                                 |  |

Rev. 2 -10- OP-EN-02

#### 1.3. 8 Electrical Safety

Please check the following conditions before attempting to use the product.

- Are you using the proper power source line? (100-240VAC)
- Are all parts (power cord and optional parts) connected properly to the product?



#### Warning

The AC power plug is a means to isolate its circuits electrically from the supply mains on all poles simultaneously. Do not place the device in an area when there is difficult to disconnect from the supply mains.



#### Warning

To avoid risk of electric this equipment must only be connected to a supply mains with protective earth.



#### Warning

The ACCURO pulse oximeter is a prescription device and is to be operated by qualified personnel only.



#### Caution

To prevent noise, install the product away from generators, X-ray machines, speakers and power cords. Proximity of the product to such equipment may result in improper functioning of the product and lead to undesirable results. A separate power circuit and secure grounding of the product are very important. Sharing of a power source with other equipment(s) may lead to undesirable results.

#### 1.3. 9 Maintenance And Cleaning

- Using various methods can clean pulse oximeter and its accessories. Please follow the methods mentioned below to avoid unnecessary damage or contamination to the Equipment.
- In the event that harmful (unauthorized) materials are used for cleaning, the damaged or contaminated Equipment shall not be serviced without charges regardless of warranty period.

Rev. 2 -11- OP-EN-02

## (Laution

Please check carefully both frame and sensor, after cleaning the Equipment, Do not use the Equipment that is worn out or damaged.

- At least once a month, clean and wipe off the frame by using the soft cloth after wetting it with lukewarm water and alcohol. Do not use lacquer, thinner, ethylene, or oxides, which could be harmful to the Equipment.
- For surface cleaning and disinfection of the monitoring system, follow institutional procedures or the recommended actions below.
- Surface cleaning Use a soft cloth dampened with either a commercial, nonabrasive cleaner or a solution of 70% alcohol in water, lightly wiping the surfaces of the monitoring system.
- Disinfection Use a soft cloth saturated with a solution of 10% chlorine bleach in tap water, lightly wiping the surface of the monitoring system.
- Before attempting to clean a Charmcare sensor, read the Instructions for Use enclosed with the sensor. Each sensor model has cleaning instructions specific to that sensor. Follow the sensor cleaning and disinfecting procedures in the particular sensor's Instructions for Use.
- Make sure both cables and accessories are free of dust or contaminants, and wipe them off with soft cloth wetted with lukewarm water (40 °C/104°F), and at least once a week, clean them by using the clinical alcohol.
- Do not submerge the accessories under any liquid or detergent. Also, make sure any liquid not to penetrate into the Equipment or probe.

## (Caution

Do not dispose single use probe to any hazard place, Always think about environmental contamination

Perform the following checks every 12 months.

- Inspect the monitoring system for mechanical and functional damage or deterioration.
- Inspect the internal fuses for proper value and rating.
- Ensure all user interface items, cables, and accessories function normally.
- At least once a month, clean and wipe off the frame by using the soft cloth after wetting it with lukewarm water and alcohol. Do not use lacquer, thinner, ethylene, or oxides, which could be harmful to the Equipment.

Rev. 2 -12- OP-EN-02



## Warning

A functional tester cannot be used to assess the accuracy of a pulse oximeter probe or a pulse oximeter monitor

#### 1.3. 10 Classifications

The ACCURO patient monitor is classified, according to IEC 60601-1 as:

| Equipment Classification      | ISO 80601-2-61 : 2012   |
|-------------------------------|---|
| Equipment Classification      | EN 60601-1 : 2012   |
| Type of protection            | Class II  |
| Degree of protection          | BF-Type applied part  |
| Electromagnetic Compatibility | EN 60601-1-2:2007   |
| Liquid Ingress                | IPX1  |
| Degree of Safety              | Not suitable for use in the presence of flammable anesthetics |

## 1.3. 11 Description Of Product & Label Symbols

| Icon               | Comments  | lcon                  | Comments  |
|--------------------|---|-----------------------|---|
| ()                 | System power On/Off                             | Ø                     | Alarm pause   |
| 1(0)               | Pulse Volume                                    | was a second          | Screen change   |
| ~                  | AC Input Indicator                              | 100-240VAC<br>50/60Hz | Connect AC power cord   |
| <u>-</u> +         | Battery Input Indicator                         | RS-232<br>10101       | Use for program upgrade   |
|                    | Type BF applied part complying with IEC 60601-1 | %SPO2                 | Oxygen Saturation   |
| ***                | Manufacturer                                    | IPX 1                 | Conforms to IEC60601-1<br>sub-clause 44.6 and<br>IEC60529 standard. |
| REF                | Reference number                                | $\sim$                | Date of manufacture   |
| SN                 | Serial number                                   | Z                     | Crossed-out wheeled bin   |
| <b>C €</b><br>0120 | CE Mark   | EC REP                | EU representative   |
| 子                  | Do not use hand hooks                           | <del>**</del> *       | Keep away from water  |
|                    | Fragile, handle with care                       | <u> </u>              | This side up  |
| $\triangle$        | Caution   | <u>^</u>              | Warning   |
|                    | Refer to instruction manual / booklet           | 0                     | Prohibition   |
|                    | Note  | LOT                   | Lot number  |
| ■<br>MENU          | Menu  | Q                     | Alarm   |

Rev. 2 -14- OP-EN-02

# Part 2

# Product Summary

#### 2.1 Intended use

The intended use of Pulse oximeter ACCURO is detecting and measurement as below parameters to provide and help doctors for making figure out the patient's vital condition.

- Exact and stable Oxygen saturation of arterial blood
- Accurate Pulse Rate

Especially, Pulse oximeter ACCURO has AGC (Auto Gain Controller) function for weak pulsatile patients.

The device is to detect and control its gain automatically in accordance with patient's pulse rate strength.

Designed for hospital, transport and home use.



Hospital use typically covers such areas as general care floors (GCFs), operating rooms, special procedure areas, intensive and critical care areas within the hospital and in hospital-type facilities. Hospital-type facilities include physician office-based facilities, sleep labs, skilled nursing facilities, surgicenters, and sub-acute centers.



Intra-hospital transport includes transport of a patient within the hospital or hospital-type facility.

- Frequently used functions
  - Connect the SpO2 probe
  - Patient mode setting
  - Alarm system (according to IEC 60601-1-8:2006)
  - Measure SpO2, Pulse rate
  - Battery (Internally powered)

Rev. 2 -15- OP-EN-02

#### ■ Essential performance

- SpO2 accuracy
- Heart rate accuracy
- Alarm limit
- Battery condition alarm

#### Patient population

a) Age: newborn, infant, adult

b) Weight: not relevantc) Health: not relevantd) Nationality: multiple

e) Patient state: patient is not user - not relevant, unless patient is agitated

#### ■ Part of the body or type of tissue applied to

a) Measurement site: finger, the back of hand, arm, toe, top of the foot.

b) Condition: the area which wish to place the probe is wiped down with alcohol.

#### User profile

#### Intended conditions of use

| Considerations   | Requirement description  |  |
|--|--|--|
| Education  | Nurse who had gotten a specialized education                                 |  |
| Education  | Physician who have the national license and can treat a patient              |  |
|  | Read and understand 'westernized Arabic' numerals when written in Arial font |  |
|  | Can perceive alarm conditions and distinguish alarm priority                 |  |
| Knowledge  | Understand means of the parameter used on the equipment                      |  |
|  | Can distinguish spo2 and Pulse rate  |  |
|  | Understands hygiene  |  |
| Language • Understand how the equipment functions in the IFU |  |  |
| understanding  | Understand abbreviations used in the IFU                                     |  |
|  | A training related to how the equipment functions                            |  |
| Evaniana   | Be enough familiar with the IFU  |  |
| Experience   | Can use spo2 probe correctly   |  |
|  | Can diagnose patient states through measured parameters                      |  |

Rev. 2 -16- OP-EN-02

#### ■ Intended conditions of use

| Considerations   |   | Requirement description                                |  |
|------------------|---|--|--|
|                  |   | Professional use, not intended for home use            |  |
|                  |   | Indoor use only  |  |
|                  |   | When it is functioning it shall keep its calibration / |  |
|                  |   | precision  |  |
|                  | General   | Use environment away from x-ray equipment, any mobile  |  |
|                  |   | device   |  |
|                  |   | Battery charging and replace for maintaining normal    |  |
|                  |   | operation  |  |
| Environment      | Conditions of   | Operation distance: within 1m                          |  |
| Environment      | visibility  | Operation distance, within 1111                        |  |
|                  | Physical  | Operating environment                                  |  |
|                  |   | Temperature: 5°C to 35°C                               |  |
|                  |   | Humidity: 30% to 85% non-condensing                    |  |
|                  |   | Storage environment                                    |  |
|                  |   | Temperature: -20°C to 70°C                             |  |
|                  |   | Humidity: 10% ~ 100%                                   |  |
|                  |   | background sound pressure level: <55db in the range of |  |
|                  |   | 100Hz – 8khz   |  |
| Frequency of use | Measure patient state continuously                      |  |  |
| Mobility         | Portable medical device to be used on a resting patient |  |  |

Rev. 2 -17- OP-EN-02

## 2.2 List of Parts

#### ■ Basic Parts

| Туре        | Quantity |
|-------------|----------|
| ACCURO Unit | 1        |

#### ■ Basic Accessories

| Туре                                   | Quantity |
|--|----------|
| SpO2 Sensor(ACCY-0A0PRB Reusable type) | 1        |
| Power cord                             | 1        |
| User Manual                            | 1        |

#### Optional Parts

| Туре |  |  |
|------|--|--|
| 1    | SpO2 Sensor(Pediatric, Neonate)        |  |
| 2    | SpO2 Sensor(Pediatric, Neonate Y-type) |  |
| 3    | SpO2 extension cable(9pin)             |  |
| 4    | IV Pole Clamp                          |  |
| 5    | Carrying Bag                           |  |

#### Sensor and Cable Length

| Туре                  | Length   |
|-----------------------|----------|
| SpO2 adult Sensor     | 1m, 3m   |
| SpO2 pediatric Sensor | 1m, 3m   |
| SpO2 neonate Sensor   | 1.5m, 3m |
| SpO2 neonate Sensor   | 1m, 3m   |
| SpO2 extension cable  | 2.5m     |
| Power Cord            | 1.8m     |

# Warning

The use of accessories, pulse oximetry sensors, and cables other than those specified may result in inaccurate readings of the monitoring system and increased emission of the monitoring system.

Rev. 2 -18- OP-EN-02

## 2.3 Nomenclature of the Parts

This section identifies the symbols, controls, displays, and indicators on the ACCURO

#### Unit

#### ■ Front Panel




| No | Name                  | Description  |  |  |
|----|-----------------------|--|--|--|
|    |                       | Connect the patient cable to the ACCURO by patient cable         |  |  |
|    |                       | plugging the cable into the patient cable connector.             |  |  |
| 1  | SpO2 probe port       | Use only connector charmcare compatible sensors and cables       |  |  |
|    |                       | with this oximeter.  |  |  |
|    |                       | See part5, sensors and patient cables for more details.          |  |  |
|    |                       | Measured SpO2 & Pulse rate values are displayed on the LCD       |  |  |
| 2  | Screen                | screen. Additional information such as waveform, alarm           |  |  |
|    |                       | messages and short trend records are displayed.                  |  |  |
| 3  | Alarm lamp            | The alarm lamp is illuminated when an alarm condition is active  |  |  |
| 3  | Alaitti lattip        | and the Alarm Status Indicator is shown.                         |  |  |
|    | Alarm pause<br>button | Press the alarm pause button to temporarily silence patient and  |  |  |
|    |                       | low battery alarms. Press the alarm pause button when the        |  |  |
|    |                       | "SENSOR OFF" or "FINGER OFF" messages are flashing (i.e.         |  |  |
|    |                       | The sensor is removed from the patient) to acknowledge the end   |  |  |
| 4  |                       | of monitoring. In these states, all further alarms are suspended |  |  |
|    |                       | until the pulse oximeter starts measuring SpO2 and pulse rate    |  |  |
|    |                       | again.   |  |  |
|    |                       | Note: System failure alarms can be paused by pressing the        |  |  |
|    |                       | power/standby or alarm pause button. If the power/standby        |  |  |

Rev. 2 -19- OP-EN-02

|          |                          | button does not pause the system fault alarm, press the alarm   |  |  |
|----------|--------------------------|---|--|--|
|          |                          | pause button.   |  |  |
| 5        | Screen Button            | This button switches from one main screen mode to another.      |  |  |
| <u> </u> | Screen bullon            | Main screen modes can be found in this user manual at Part 6.   |  |  |
| 6        |                          | Used to move upward from one menu category to another or to     |  |  |
| 0        | <b>∧</b> button (up)     | alter available preset values of each menu category.            |  |  |
|          |                          | Change the screen contrast ① at the top menu                    |  |  |
| 7        | > button (right)         | Used to move rightwards from one menu category to another or    |  |  |
|          |                          | to alter available preset values of each menu category.         |  |  |
| 8        | M. hutton/down           | Used to move downward from one menu category to another or      |  |  |
|          | V button(down)           | to alter available preset values of each menu category.         |  |  |
|          | <b>&lt;</b> button(left) | Change the pulse volume (a) at the top menu                     |  |  |
| 9        |                          | Used to move leftwards from one menu category to another or to  |  |  |
|          |                          | alter available preset values of each menu category.            |  |  |
| 10       | selection                | Used to move up and down on the main menu tree and to select    |  |  |
|          | button                   | each menu category.   |  |  |
| 11       | AC power usage           | Indicates connection to alternating current power source.       |  |  |
| ''       | indicator                | indicates connection to alternating current power source.       |  |  |
| 12       | Battery status           | The Battery Status Indicators show the capacity.                |  |  |
| 12       | indicator                | The Dattery Status indicators show the capacity.                |  |  |
|          |                          | Press the Power On/Off Button to turn the instrument on.        |  |  |
| 13       | Power button             | Press, hold the button for more than 2 seconds and then release |  |  |
|          |                          | the button to turn the instrument off.                          |  |  |

Rev. 2 -20- OP-EN-02

#### Rear Panel



| No | Name                  | Description   |  |  |
|----|-----------------------|---|--|--|
|    |                       | It can directly supply AC power via this AC power port for      |  |  |
| 1  | AC Power port         | supplying constant power and recharging the handheld while      |  |  |
|    |                       | docked.   |  |  |
| 2  | Serial port           | This ports are software upgrade port                            |  |  |
|    | Connector             |   |  |  |
| 3  | Attaching clamp holes | Holes for attachment to the IV pole clamp                       |  |  |
| 4  | Speaker               | The speaker indicates audio alarms. Care should be taken not to |  |  |
| 4  |                       | cover the speaker and muffle the audible alarm volume.          |  |  |



## Warning

To prevent electric shock, keep the top cover closed and do not attempt to disassemble the product on your own. Product should only be disabled by authorized service personnel.

Rev. 2 -21- OP-EN-02

## Part 3

## **Product Installation**

First-time users: Read the installation instructions thoroughly and install the product in a safe place to ensure product longevity.

#### **Caution When Installing**

- Use the product in a place where the temperature is between  $5\sim35$  °C and the humidity is  $30\sim85$ %.
- Make sure that the power cord is plugged in properly.
- Do not plug in more than one device into one outlet.
- Ground the product if noise develops.
- Do not use a power cord that causes noise.
- The product is very sensitive to shock. Caution is advised.
- Keep the product dust-free and install away from flammable materials.

#### **Connecting The AC Power**

- Connecting the Power Connect the power to the power port.
- Grounding the Product Connect the grounding line to the grounding port in the back of the product if needed.

AC power cord is connected in AC power port, a green light appears in the AC power usage indicator. Product usage is possible without AC power via the embedded battery. In such cases, a red light appears in the Battery Usage Indicator.













② Battery Usage Indicator

## **A** Caution

When the charge level of the embedded battery is running low, the battery usage indicator blinks and an information alarm is emitted. Depending on usage parameters, power may be cut instantly, and thus the AC power should be connected for use.

Rev. 2 -22- OP-EN-02



#### Caution

Both the monitoring system and the instrument connected to it must utilize a grounded outlet.



#### Warning

The battery usage warranty period is six months. Following this period, battery life can diminish significantly.



#### Warning

Use the power supply device supplied by the manufacturer. Failure to do so may result in electric shock and damage to the product. Edge or points of the Product can hurt to patient and user.




#### Warning

Accuracy may degrade if it is connected to secondary I/O devices when the equipment is not connected to earth reference.

#### Installing the IV Pole


The product may also be mounted on a clamp, which can then be attached to an IV Pole. This makes the product more portable. (\*IV Pole is optional.)

(1) The product may be used horizontally or vertically. First, decide how you want the product to be mounted. Then, insert two bolts into the holes on the clamp and tighten.



Holes for attachment to the stand

(2) Mount the clamp to the pole by sliding the clamp onto the pole and then, turn the handle to tighten.





The pole on which the stand is mounted should be less than 30mm in diameter. Do not use poles that are any thicker than this.

Rev. 2 -24-OP-EN-02

# Part 4

# **Using The Product**

This section contains information on the basic nomenclature and directions for using the product.


#### **Preparations before Use**

- Check for defects on the exterior of the SpO2 sensor(Oxygen sensor).
- Check for defects on the exterior of the AC power cord, oximetry use cable.
- Connect the power, oximeter use cable to the SpO2 sensor and product. (If the oximetry use cable is not utilized, SpO2 sensor is directly connected to the product.)
- Check the connection status of the cable and status of SpO2 sensor before turning the product on and connecting the product to the patient.
- Since this product uses the principle of a spectrophotometer, check if a surgery light is located where the product is placed, and if the ray blocks the usability cover the SpO2 sensor with a opaque cloth.
- Check the attachment region of the SpO2 sensor and the compatibility, and avoid certain locations according to the status and use a SpO2 sensor that can be attached on to the bridge of the nose or forehead.
- When attaching the SpO2 sensor on to the finger, no manicure should be on the finger nails.
- For sterilization or antiseptic. cleaning is necessary for the SpO2, select the methods shown below according to the environment conditions of the medical institution equipment, and use the usage and application of the professional prefers.

<sup>\*</sup>How to use antiseptic .cleaning solution

| Chemical Product |                   | Method                  |
|------------------|-------------------|-------------------------|
| Alcohol          | Isopropyl alcohol | Use soaked gauze.       |
| Liquid Soon      | Benzalkonium      | 0.05 W/V%(x200 diluted) |
| Liquid Soap      | Chloride          | Use soaked gauze.       |
| lodine           | Povidone-iodine   | 0.02 W/V%(x50 diluted)  |
|                  |                   | Use soaked gauze.       |
| Glutaral         | C5H8O2            | 2 W/V %(test solution)  |
| Giularai         | G31 1002          | Use soaked gauze.       |


Rev. 2 -25- OP-EN-02



| No | Name                 | Description   |
|----|----------------------|---|
| 1  | SpO2/PR data area    | Displays the measured SpO2, PR and related values.            |
| 2  | Waveform area        | Displays SpO2 waveform.                                       |
|    |                      | Displays measured perfusion Index.                            |
| 3  | PI data area         | The Perfusion Index indicates numerically the percentage of   |
|    |                      | pulsatile signal to non-pulsatile signal (Pulse strength).    |
| 4  | Patient Alarm area   | Displays the alarm message from the measured value. When      |
|    | T attent Alaim area  | multiple alarms occur, it shows according to higher priority. |
| 5  | Technical Alarm area | Displays the alarm message related to the system.             |
| 6  | Short Trend area     | Displays the SpO2, PR value saved during the previous 30      |
| U  | Siloit Hellu alea    | minutes or 15 minutes on a graph.                             |
| 7  | System status area   | Displays status of battery, alarm and current time.           |

Rev. 2 -26- OP-EN-02

## SpO2/PR Data Area



| No  | Name              | Description   |
|-----|-------------------|---|
| 1   | Oxygen saturation | Displays the measured value of SpO2. When the value exceeds the     |
| '   | (SpO2) reading    | SpO2 alarm limit, the SpO2 value will turn red and flash.           |
|     |                   | The Saturation Alarm Limits Display shows the upper and lower       |
| 2,3 | SpO2 alarm limit  | saturation alarm limits. When the value above or below limits if it |
|     |                   | turns red.  |
| 4   | SpO2 alarm status | Displays the alarm watch state of SpO2. (refer to alarm status)     |
| 5   | Pulse rate (PR)   | Displays the measured PR value. When the value exceeds the PR       |
| 5   | reading           | alarm limit, the PR value will turn red and flash.                  |
|     |                   | The pulse rate Alarm Limits Display shows the upper and lower       |
| 6,7 | PR alarm limit    | saturation alarm limits. When the value above or below limits if it |
|     |                   | turns red.  |
| 8   | PR alarm status   | Displays the alarm watch state of PR. (refer to alarm status)       |
| 9   | Heart beat icon   | Blinks when a heart beat occurs. If the PR value derails from PR    |
| 9   |                   | alarm limit it will be blinking in red.                             |
|     | Level bar         | Indicates pulse beat and the relative (non-normalized) pulse        |
| 10  |                   | amplitude in numbers only view.                                     |
| 10  |                   | As the detected pulse becomes stronger, more bars light with each   |
|     |                   | pulse.  |

Rev. 2 -27- OP-EN-02

#### **System Status**



| NO | Name               | Description  |  |  |
|----|--------------------|--|--|--|
|    |                    | Displays the battery status.                             |  |  |
|    | Power status       | Uses a battery and capacity is low(about10min)           |  |  |
|    |                    | Uses a battery and capacity is low (about30min)          |  |  |
|    |                    | Uses a battery and capacity is 1/4                       |  |  |
| 1  |                    | Uses a battery and capacity is 1/2                       |  |  |
|    |                    | Uses a battery and capacity is 3/4                       |  |  |
|    |                    | Uses a battery and capacity is FULL                      |  |  |
|    |                    | Uses AC power and battery is being charged               |  |  |
|    |                    | Uses AC power only but battery is fully charged          |  |  |
|    |                    | Displays the alarm sound off state.                      |  |  |
| 2  | Alarm sound status | 0:48 Displays the alarm sound paused state and left time |  |  |
|    |                    | 🔉 0 Displays the Alarm sound off state                   |  |  |
| 3  | Pulse volume       | Display pulse volume level and state                     |  |  |
| 4  | LCD brightness     | Display brightness level                                 |  |  |
| 5  | Patient type       | Display patient type  ♣ ADULT  ♣ PEDIATRIC  ♣ NEONATE    |  |  |
| 8  | Current-time clock | Displays the current time(hour: minute)                  |  |  |

## **Post-Use Storage And Management**

- After use, store the product in an area that does not exceed 85% humidity, and maintain the temperature in the range of 5°C~ 35°C.
- Do not store the product together with chemical products and/or in an area subject to gas exposure.
- Avoid exposing the product directly or indirectly to heat.
- Make sure to review the product usage instructions prior to use.
- Store the product in an area where it will not come in contact with water.
- Store the product in an area free of dust and/or other foreign substances.
- Store the product away from direct sunlight.

Rev. 2 -28- OP-EN-02

- Maintain the product in a safe and stable condition away from vibration and/or other sources of shock.
- If the product and/or sensor area becomes dirty, wipe it clean with rubbing alcohol and cotton swabs, and then allow it to dry at room temperature for 30 minutes.

#### Power off

 Press power button during 3seconds to power off. When display is off, separate the power cord by pulling off the cords connector from device rear.


Rev. 2 -29- OP-EN-02

## Part 5

# Measuring the SpO2

#### **Attaching The SpO2 Probe**

- (1) Wipe down the area where you wish to place the probe with alcohol.
- (2) Attach the probe to the patient's finger.
- (3) To get an accurate reading, make sure the patient minimizes all movement and Please attach the probe wire to patient's finger firmly. Attach the bandage loosely so as to not cut of circulation to the finger.
- (4) Check the patient's finger and the probe every two to three hours to make sure that the sensor is properly placed over the finger. If there is a change in the appearance of the patient's finger due to the prolonged exposure to the probe, switch the probe to another finger.





#### Warning


- -Measuring SpO2 on a patient undergoing an MRI may result in severe burns for the patient. To minimize risk for burns, use a non-inductive wire. In the event that this does occur, immediately remove the probe from the patient.
- -The area around the Sp02 sensor shall not exceed 37  $^{\circ}$ C . The sensor will not work in temperatures above 37  $^{\circ}$ C .
- -Do not attach the probe near arterial or venous catheters.
- -Make sure that the sensor emits a light and that the sensor is properly placed over the patient's finger.
- -Excessive pressure for prolonged periods from the sensor may cause necrosis of the skin.
- -Sensor sites should be changed at least every 8 hours

Rev. 2 -30- OP-EN-02

#### **SpO2 Sensor Port**

Use only approved sensors and interface cables when connecting to the sensor port. Connecting any other cable or sensor influences the accuracy of sensor data, which may lead to adverse results. If the pulse beep tone does not sound with each pulse, the pulse beep volume is set to zero, the speaker is malfunctioning, or the signal is corrupt. Reset the device.

#### - Wire connection





-Before use the probe sensor, operator needs to verify that the SpO2 sensor type match

## (Caution

- -Handle the probe sensor and wire with caution. Careless handling may damage the sensitive sensor. Protect the wire from sharp objects.
- -The skin of patients who have high fevers or have problems with distal circulation will be 2-3 degrees higher than normal.
- -Patients with abnormally high oxyhemoglobin or methemoglobin levels will not give a proper SpO2 reading.

Rev. 2 -31- OP-EN-02



- -Taking the NIBP can affect the SpO2 reading. When taking NIBP, place the SpO2 probe on the other arm.
- -Avoid using the probe with other medical equipment that affects blood flow. Avoid placing the probe near an area that requires medical attention.

#### **Measuring The SpO2**

- (1) Connect AC power to AC power port on the back side.
- (2) Connect SpO2 sensor to upper side.
- (3) Press the power button.
- (4) After initialed ACCURO, the measurements are ready.
- (5) Attach the SpO2 sensor to finger.
- (6) When the display reads "Learning," this means that SpO2 measurement has begun. Minimize patient movement until the screen no longer reads "Learning" to ensure a proper SpO2 reading.
- (7) After receiving data for 3 seconds, SpO2 results will be indicated.

Rev. 2 -32- OP-EN-02

# Part 6

## Switch Screen Mode

#### **Switch Screen Mode**

ACCURO consists of 4 types of screens and the can be switch to the screen the user designates. Press button to switch the screen. Switching the screen is as follows.




SpO2, PR, Wave, Short trend <screen 1>



SpO2, PR, Short trend <screen 2>



SpO2, PR, Wave



Wave, SpO2, PR <screen 4>


Rev. 2 -33- OP-EN-02

# Part 7

## How to Use Menu

The SETUP menu allows for the adjustment of various settings including current time, patient mode, alarm volume, pulse volume, alarm settings, etc.

#### **Entering Main Menu**



- (1) Press the menu button to enter the SETUP menu. Screen will change to the right of the <screen 4> to display the parameter data is the right of the menu screen.
- (2) Move the icon to enter the sub-menu <(left) button or > (right) button, press the menu button.
- (3) Press the menu button on the "EXIT" item will leave the SETUP menu.

#### **Setup Configuration**



- (1) Press the  $\Lambda$  (up) button or V (down) button to move to the item set.
- (2) You can adjust the settings in the item set press < (left) button or > (right) button.
- (3) Press the menu button on the "EXIT" item will leave the SETUP menu.

Rev. 2 -34- OP-EN-02

The contents of each item are as shown as below.

| Item          | Description                                    | Range             |
|---------------|--|-------------------|
| PULSE VOLUME  | Set the level of the pulse volume.             | OFF, 1~7          |
| BUTTON VOLUME | Set the level of the button volume.            | OFF, 1~7          |
| BRIGHTNESS    | Set the brightness of the screen.              | 1~5               |
| SCREEN TYPE   | Set the main screen mode.                      | 1~4               |
| SPO2 AVERAGE  | Set the calculation time for SpO2 measurement. | 4, 8, 12, 2 (sec) |
| SPOZ AVERAGE  | Set the calculation time for SpO2 measurement. | (Neonate 8,12)    |
| WAVE FILL     | Set whether to fill in the wave at display.    | OFF, ON           |
| PLETH BAR     | Set the display status of the level bar.       | OFF, ON           |
| SWEEP SPEED   | Set display speed for ECG wave.                | 6.25, 12.5,       |
|               |  | 25mm/s            |
| TIME          | Go to submenu to set time                      |                   |

## **Setup Current Time**



- (1) From the setup menu using menu buttons move to time submenu.
- (2) Press the  $\wedge$  (up) button or  $\vee$  (down) button to move to the item set.
- (3) You can adjust the settings in the item set press < (left) button or > (right) button.
- (4) Press the menu button on the "EXIT" item will leave the SETUP menu.

The contents of each item are as shown as below.

| Set   | Year      | Month | Day  | Hour | Minute |
|-------|-----------|-------|------|------|--------|
| Range | 2000~2099 | 1~12  | 1~31 | 0~23 | 0~59   |

Rev. 2 -35- OP-EN-02

#### **Setup Alarm**



- (1) Press the  $\Lambda$  (up) button or V (down) button to move to the item set.
- (2) You can adjust the settings in the item set press < (left) button or > (right) button.
- (3) Press the menu button on the "EXIT" item will leave the SETUP menu.


The contents of each item are as shown as below.

| Item             | Description   | Range       |  |
|------------------|---|-------------|--|
|                  | Set the level of the alarm volume. If you want to set the |             |  |
| ALARM VOLUME     | alarm volume at "0" you can only do this by executing     | (OFF), 1~7  |  |
|                  | "Alarm 0 USE" command at the MAINT menu.                  |             |  |
| ALARM DURATION   | Set the duration (second) when in alarm sound pause and   | 30, 60, 90, |  |
| ALARIVI DURATION | alarm sound   | 120 seconds |  |
|                  | Checks if the SPO2 measurement value is out of the        |             |  |
| SPO2 ALARM       | range of the setup range during alarm watch. SpO2 alarm   | OFF, ON     |  |
|                  | watch status displays 💥 when setting is 'OFF'.            |             |  |
|                  | Checks if the PR measurement value is out of the range    |             |  |
| PR ALARM         | of the setup range during alarm watch. PR alarm watch     | OFF, ON     |  |
|                  | status displays 💥 when setting is 'OFF'.                  |             |  |
| SpO2 HIGH        | Sets the alarm upper limit for the SPO2 value. When the   | 20~100 OEE  |  |
| эрог півп        | value is larger or the same, the alarm beeps.             | 20~100, OFF |  |
| C=03   OW        | Sets the alarm lower limit for the SPO2 value. When the   | OFF 20-100  |  |
| SpO2 LOW         | value is smaller or the same, the alarm beeps.            | OFF, 20~100 |  |
| DD LIICH         | Sets the alarm upper limit for the PR value. When the     | 15~300, OFF |  |
| PR HIGH          | value is larger or the same, the alarm beeps.             |             |  |
| DD LOW           | Sets the alarm upper limit for the PR value. When the     | OFF 15, 200 |  |
| PR LOW           | value is smaller or the same, the alarm beeps.            | OFF, 15~300 |  |

# (Caution

Cautiously setup value when the user cannot recognize the patient status when the alarm occurs and alarm sound is OFF. Frequently check the patient status when OFF.

#### **Setup Patient**



- (1) Press the  $\wedge$  (up) button or  $\vee$  (down) button to move to the item set.
- (2) You can adjust the settings in the item set press < (left) button or > (right) button.
- (3) Press the menu button on the "EXIT" item will leave the SETUP menu.

The contents of each item are as shown as below.

| Item         | Description               | Range                   |
|--------------|---------------------------|-------------------------|
| PATIENT TYPE | Sets up the patient mode. | ADULT,PEDIATRIC,NEONATE |

When the Patient Mode changes, the alarm limit value and SPO2 average switches to the USER CONFIG value or FACTORY SET value.

#### **Entering Maintenance Menu**



- (1) From the SYSTEM menu using <, > arrow buttons move to MAINT, press menu button MAINT PASSWORD submenu to pop up as above.
- (2) You can press the Menu button to go to each place and change the gray background to change the number. Press the Menu button again to enter Password 4 of the value, you move from the '<' button to move to "MAINT" menu.

Rev. 2 -37- OP-EN-02

- (3) Incorrect PASSWORD is displayed at the bottom of the input "Please enter the correct password" message.
- (4) Press the menu button on the "EXIT" item will leave the SETUP menu.



MAINT setting can be manipulated by authorized personnel only. Be aware that the contents of MAINTENANCE are of nature that is not to be altered by normal users.

#### **Setup Maintenance**



- (1) Press the  $\wedge$  (up) button or  $\vee$  (down) button to move to the item set.
- (2) You can adjust the settings in the item set press < (left) button or > (right) button.
- (3) Press the menu button on the factory set load items displayed "factory set loaded" message and press the menu button on the user config save items displayed "user configuration saved." message.
- (4) Press the menu button on the "EXIT" item will leave the SETUP menu.

The contents of each item are as shown as below.

| Item             | Description   | Range   |  |
|------------------|---|---------|--|
| ALARM 0 USE      | From SOUND menu it enables to turn ALARM Volume 0.        | ON, OFF |  |
| ALARM 0 REMIND   | When the alarm volume is 0 and the alarm occurs, sets if  | ON, OFF |  |
| ALAKINIO KLININD | the alarm should beep in 3 minutes duration.              | ON, OFF |  |
| LANGUAGE         | Changes language setting, as altered once with newly      | ENGLISH |  |
| LANGUAGE         | desired language it will display it immediately.          | ENGLISH |  |
| FACTORY SET      | Press the Menu button to change the factory default       | LOAD    |  |
| FACTORT SET      | settings, all settings.                                   | LOAD    |  |
| USER CONFIG      | The user-defined value is stored in memory, brings to the | SAVE    |  |
| USER CONFIG      | default settings when you reboot.                         | SAVE    |  |

Rev. 2 -38- OP-EN-02

### Part 8

### Save and View TREND Data

#### **Saving Measured Data**

The value of the measured SpO2, PR value is saved as below.

- 1. The SpO2, PR value is saved in a 2, 4, or 10 second unit span.
- 2. Saves the limit alarm state.
- 3. When data exceeds the capacity, old data is deleted and new data is saved.
- 4. When the SpO2 sensor is not connected to the device, data is not saved.

#### **Setup Trend**



- (1) Press the  $\Lambda$  (up) button or V (down) button to move to the item set.
- (2) You can adjust the settings in the item set press < (left) button or > (right) button.
- (3) Press the menu button on the "EXIT" item will leave the SETUP menu.


The contents of each item are as shown as below.

| Item          | Description  | Range         |
|---------------|--|---------------|
| SAVE INTERVAL | Sets the interval of saving time.                          | 2, 4, 10 sec. |
| SHORT TREND   | Sets total duration of short trend until present time from | 15, 30 min    |
| SHORT TREND   | the set value ago. (For some available main screens)       |               |
| ERASE DATA    | Erase trend data.  |               |


Rev. 2 -39- OP-EN-02

#### **Short Trend View**

The recently recorded SpO2 and PR value can be shown in a graph as shown below on the Available on main screen type 1 and 2 shows SHORT TREND up till present time from the preset value (SETUP>TREND>SETUP>SHORT TREND)



#### **Tabular Trend View**



- (1) Press menu button on the SETUP>TREND>TABULAR.
- (2) TABULAR trend table pops up like below. It show the most recent data.

| Item             | Description   |
|------------------|---|
| ∧ (up) button    | Go to the previous page. Moves quickly to the page after pressing more than 2 seconds   |
| V (down) button  | Go to the next page. Moves quickly to the page after pressing more than 2 seconds < (left) button Displays the page of the newest saved data. |
| < (left) button  | Go to the first page.   |
| > (right) button | Go to the last page.  |
| MENU button      | Exit tabular trend menu. Moves to the upper menu (TREND) screen.  |


When saving a pair of parameter values (SpO2 and PR), it saves any alarm events as well. (If any alarm ever goes off at any given saving time – display red color)

Rev. 2 -40- OP-EN-02

### **Graphic Trend View**



- (1) Using <, > arrow button move to items (30min~72hr) and press MENU button.
- (2) GRAPHIC trend pops up like below.



| No  | Item            | Description  |
|-----|-----------------|--|
| 1   | Time interval   | Displays the size of time of the Graphic trend screen.   |
| 2   | Current page    | Moves to the latest data saved.  |
| 3,4 | Move pages      | Using ARROW(A, V) buttons you can move along pages.  |
| 5   | Move trend bars | Using ARROW(<, >) buttons you can move along trend bars.   |
| 6   | Trend bar       | Displays current trend bare on present page.   |
| 7   | Date/time       | Date and time of trend data saved on the timeline.   |
| 8   | SpO2 value      | Displays the measured SpO2 value in the trend bar. The value color changes to red when the limit alarm occurs. |
| 9   | PR value        | Displays the measured PR value in the trend bar. The value color changes to red when the limit alarm occurs.   |
| 10  | EXIT            | Moves to the upper menu (TREND) screen using the 'MENU' button.  |

Rev. 2 -41- OP-EN-02

### Part 9

### Alarm

The screen, sound and lamp alerts according to the alarm sensed alarm signaled from the patient during measurement and alarm occurred from the system

#### **Alarm Categories**

Alarm consists and is managed in two types, Patient alarm and Technical alarm.

Patient Alarm

This alarm occurs from the patient data. Occurs when the measured physiological signal value initiate the limit alarm or alerts through the patient status.

Technical Alarm

This alarm occurs from the system. Occurs for product related actions or system errors.

General Message

The alarm does not occur and shows a related message to the user. This is information such as "LEARNING..." and "PATIENT MOVING" that starts the measurement after sensing a signal.

#### Alarm Levels

Single alarm is classified and managed in high, medium and low priority level.

High Priority

Occurs when the danger is detected in life of the patient.

Medium Priority

Occurs when the physiological signal is abnormal or requires treatment.

Low Priority

Occurs when the user has to be notified or the system confronts an error during measurement.

Patient alarm and Technical alarm are prioritized according to the content, and some parts can be set by the user from alarm settings.

Note . Refer to Attachment2 for alarm message and alarm level content.

#### **Alarm Indicators**

The alarm is displayed as below when alarm occurs.

Rev. 2 -42- OP-EN-02

#### Alarm message

In the upper portion of the screen the alarm message field shows as below.



Patient alarm message area
 Displays the Patient alarm message.
 Technical alarm message area
 Displays the Technical alarm message.

Displays the message according to high priority, and shows the message with same priority in a sequence.

High Priority: Red color font, and displays patient alarm attaches \*\*\*.

Medium Priority: Yellow color font, and displays patient alarm attaches \*\*.

Low Priority: Yellow color font, and displays patient alarm attaches \*.

#### Alarm Tone

When multiple alarms occur, alarms will beep according to the order of top priority.

When the state is on alarm sound pause or alarm sound off, alarm sound will not occur.

| Tone Category                  | Description                          |
|--------------------------------|--------------------------------------|
| ***High priority alarm signal  |                                      |
| Pulse Frequency                | 523Hz                                |
| Pulse width                    | 150ms                                |
| Number of pulse in burst       | 10, interburst interval of 4 seconds |
| Repetitions                    | continually                          |
| Sound level range              | 70~76dB                              |
| **Medium priority alarm signal |                                      |
| Pulse Frequency                | 523Hz                                |
| Pulse width                    | 150ms                                |
| Number of pulse in burst       | 3, interburst interval of 8 seconds  |
| Repetitions                    | continually                          |
| Sound level range              | 60~67dB                              |
| *Low priority alarm signal     |                                      |
| Pulse Frequency                | 659Hz                                |
| Pulse width                    | 250ms                                |

Rev. 2 -43- OP-EN-02

| Number of pulse in burst | 2, interburst interval of 16 seconds |
|--------------------------|--------------------------------------|
| Repetitions              | continually                          |
| Sound level range        | <b>57~62</b> dB                      |

#### Alarm lamp


When multiple alarms occur, lamps will turn on according to the order of top priority.

High Priority : Red lamp blinks 2 times for 1 second.Medium Priority : Yellow lamp blinks for each second.


Low Priority : Yellow lamp is turned on and maintains its state.

#### **Priority Signal Timing**


#### 1) High priority



#### 2) Medium priority



#### 3) Low priority



Rev. 2 -44- OP-EN-02

#### **Alarm Conditions**

| Alarm condition                   | Auditory alarm | Alarm signal    | Messages       |  |
|-----------------------------------|----------------|-----------------|----------------|--|
| High heart rate limits violated   | Medium         | Visual alarm    | HIGH PR        |  |
| riigii neartrate iimits violateu  | Mediairi       | Auditory alarm  | HIGHTIK        |  |
| Low heart rate limits violated    | Medium         | Visual alarm    | LOW PR         |  |
| Low Healt Fate IIIIII Worlded     | Wicalani       | Auditory alarm  | LOWIN          |  |
| High SpO2 limits violated         | Medium         | Visual alarm    | HIGH SPO2      |  |
| riigii opoz iiriito violatea      | Wicararr       | Auditory alarm  | 1110110102     |  |
| Low SpO2 limits violated          | Medium         | Visual alarm    | LOW SPO2       |  |
| Low opez innits violated          | Wicalani       | Auditory alarm  | LOW 01 02      |  |
| Unable to determine pulse rate or | Low            | Visual alarm    | CHECK PROBE    |  |
| SpO2.                             | LOW            | Auditory alarm  | ONEOKT KOBE    |  |
| Communication Error               | High           | Visual alarm    | SPO2 FAULT     |  |
| Communication Error               |                | Auditory alarm  |                |  |
| SpO2 cable/sensor disconnect      | Low            | Visual alarm    | SENSOR OFF     |  |
| Opoz dabie/seriodi disconinect    |                | Auditory alarm  |                |  |
| Finger off form SpO2 sensor       | Low            | Visual alarm    | FINGER OFF     |  |
| r inger on form opez senser       |                | Auditory alarm  |                |  |
| Low battery                       | Medium         | Visual alarm    | LOW BATTERY    |  |
| Low Ballery                       |                | Auditory alarm  | LOW BATTERT    |  |
| Low Internal Coin Battery         | Low            | Visual alarm    | BATT ERROR     |  |
| Low Internal Com Battery          |                | Auditory alarm  | B/(IT EI(I(O)) |  |
| Clock settings lost               |                | Visual alarm    |                |  |
| (Invalid date and time value      | Low            | Auditory alarm  | CLOCK ERROR    |  |
| detected during start up)         |                | , aditory diami |                |  |

#### **Alarm Conditions Delay**

#### Average Time

Average time 2 - Responds to changes in blood oxygen saturation in 2 to 4 seconds.

Average time 4 - Responds to changes in blood oxygen saturation in 4 to 8 seconds.

Average time 8 - Responds to changes in blood oxygen saturation in 8 to 16 seconds.

Average time 12 - Responds to changes in blood oxygen saturation in 12 to 24 seconds.

Rev. 2 -45- OP-EN-02

#### **Alarm Status**

#### Alarm Status Symbol

Displays the alarm watch state as shown below.

| No. | Status             |
|-----|--------------------|
| Q   | Alarm watch on     |
| Φ   | Alarm sound paused |
| X   | Alarm sound off    |
| ×   | Alarm watch off    |

#### Alarm sound paused

Press the button when the alarm occurs and the alarm sound is paused. The alarm does not sound when it is paused. During the audible alarm pause period, you can press Alarm button again to re-enable the audible alarm tone. Alarm pause time is set to SETUP > ALARM > ALARM DURATION.

There are operations according to the alarm when alarm sound is paused.

■ Patient Alarm : Stops the alarm sound during ALARM DURATION

■ Technical Alarm: When SENSOR OFF, FINGER OFF, LOW BATTERY occurs,

Alarm sound stops and is maintained before each state is canceled.



#### Warning

Do not pause button the audible alarm or decrease its volume if patient safety could be compromised.



#### Warning

Each time the monitor is used, check alarm limits to make sure that they are appropriate for the patient being monitored.

#### Alarm off

This turns the SpO2 or PR limit alarm watch off. It is excluded in alarm messages and alarm sound operations for limit alarms.

To stop the limit alarm watch of SpO2, uncheck the SETUP > ALARM > SPO2 ALARM. And the alarm watch status switches to **X** in the SpO2 screen window.

To stop the limit alarm watch of PR, uncheck the SETUP > ALARM > PR ALARM. And the alarm watch status switches to  $\bowtie$  in the PR screen window

Rev. 2 -46- OP-EN-02

#### Alarm sound off

This turns the alarm souund off. It operates the alarm message and alarm lamp when an alarm occurs.

When the value is set to '0' in SETUP > ALARM > ALARM VOLUME, the alarm souund stops, and the alarm watch status of SpO2 and PR switches to  $\bigcirc$  0.

When the alarm souund is off, an alarm sounds every 3 minutes to alert the current alarm status. This can be set when the alarm occurs and when SETUP > MAINT > ALARM 0 REMIND is checked.



#### Warning

Ensure the speaker is clear of any obstruction. Failure to do so could result in an inaudible alarm tone.



#### Warning

A potential hazard can exist if different alarm presets are used for the same or similar equipment in any single area



#### Warning

Auditory alarm signal sound pressure levels that are less than ambient levels can impede operator recognition of alarm conditions



#### Note

As the alarm sounds emitted are difficult to describe in words, it is recommended that the user listen directly to the sounds of the high priority and medium priority alarms prior to product use



#### Caution

When turning alarm sound OFF, you should take extra caution because there is a chance that you can not recognize the status of patient. If you have to turn it OFF, check on patient frequently.



#### Caution

If alarm limits set to extreme values, that can render the alarm system useless.

Rev. 2 -47- OP-EN-02

#### **Verifying Visual and Audible Alarm Indication**

If the monitor fails to perform as specified in this test, contact qualified service personnel or your local supplier for assistance.

You can verify the alarm operation for all parameters SpO2 by following the below procedures.

- (1) Connect the monitor to an AC power source.
- (2) Press Power button to turn on the monitor.
- (3) Connect the simulator to sensor input cable and connect cable to monitor.
- (4) Set the simulator to smaller value than the lower alarm limit on the monitor.
- (5) Verify following the monitor reaction:
  - a. The monitor begins to track the physiological signal from the simulator.
  - b. After about 10 to 20 seconds, the monitor displays the value measured as specified by simulator. Verify values are within the tolerances specified in Specification section for each parameter.
  - c. Audible alarm sounds.
  - d. The numerical area flashes, indicating the parameter has violated default alarm limits.



It is recommended that you listen to the high and mid alarms beforehand in order to be cognizant of the sounds when the alarms do go off.

Rev. 2 -48- OP-EN-02

# Part 10 Basic Troubleshooting

Here are some basic troubleshooting hints. There may be times when the product does not work as it is supposed to and you may not know how to proceed. Try the following steps to work through the problems.

#### General

| Situation                                  | Cause   |
|--|---|
| The power will not come on.                | - Make sure the product is plugged in.                |
|  | - Make sure the battery is fully charged.             |
| I can't hear the alarm or the pulse sounds | - Make sure the alarm is not turned off.              |
| Coin battery image (⊕) is shown in the     | - Make sure the coin battery is not lacking capacity. |
| lower portion of the screen.               | - Wake sure the com battery is not lacking capacity.  |
| The saved value is initialized.            | - Make sure the coin battery is not lacking capacity. |
|  | - Make sure the memory has no problems                |

#### Measuring the SpO2

| Situation                              | Cause                                  |
|--|--|
| The red sensor light does not come on. | - The probe is not connected properly. |
|  | - The probe is broken.                 |
| I can't see the waves.                 | - The probe is not connected properly. |

#### ■ Error Message

| Message     | Cause  |
|-------------|--|
| HIGH PR     | SpO2 is above the alarm limit. Check patient immediately.                        |
| LOW PR      | SpO2 is below the alarm limit. Check patient immediately.                        |
| HIGH SPO2   | SpO2 is above the alarm limit. Check patient immediately.                        |
| LOW SPO2    | SpO2 is below the alarm limit. Check patient immediately.                        |
| CHECK PROBE | Sensor not attached to patient. Reposition or replace sensor.                    |
| SPO2 FAULT  | Analog board resulted in unstable contact with the IC Communication. Return      |
| SPO2 FAULT  | to a qualified service technician.   |
| SENSOR OFF  | Sensor not attached to patient. Reposition or replace sensor.                    |
| FINGER OFF  | Sensor not attached to patient finger. Reposition or replace sensor.             |
| LOW BATTERY | Connect to AC power. Continue charging on AC                                     |
| BATT ERROR  | Low Internal Coin Battery. Return to a qualified service technician.             |
| CLOCK ERROR | Date and time setting lost. Set the date and time. Return to a qualified service |
| CLOCK ERROR | technician.  |

Rev. 2 OP-EN-02 -49-

## Part 11

# **Product Specification**

### **Product Specification**

#### PERFORMANCE

| Saturation (SpO2) |                  |              |             |
|-------------------|------------------|--------------|-------------|
| Range             | 0~100%           |              |             |
| Resolution        | 1%               |              |             |
| Accuracy          | Adult, Pediatric | 70%~100%     | ±2 Digits   |
|                   |                  | 50%~69%      | ±3 Digits   |
|                   |                  | 0%~49%       | unspecified |
|                   | Neonate          | 70%~100%,    | ±3 Digits   |
|                   |                  | 50%~69%      | ±4 Digits   |
|                   |                  | 0%~49%       | unspecified |
| Average           | Adult, Pediatric | 2, 4, 8, 12, |             |
|                   | Neonate          | 8, 12        |             |
| Pulse rate        |                  |              |             |
| Range             | 30~250bpm        |              |             |
| Resolution        | 1bpm             |              |             |
| Accuracy          | 30~250bpm        | ±3 Digits    |             |
| Perfusion index   |                  |              |             |
| Range             | 0.05 ~ 20%       |              |             |

#### Electrical

| Instrument       |   |  |
|------------------|---|--|
| Power Input      | 100-240 Vac, 50/60Hz                              |  |
| Fuse rating      | 5A, 32VAC/DC                                      |  |
| Battery          |   |  |
| Туре             | 7.2V Li-ion rechargeable battery                  |  |
| Charge time      | 6 hours   |  |
| Battery Capacity | Typically 6hours using new, fully charged battery |  |
| Warranty         | 6 months  |  |

Rev. 2 -50- OP-EN-02

#### Environmental

| Storage Temperature | -20℃ ~ 70℃                     |
|---------------------|--------------------------------|
| Storage Humidity    | 10% ~ 100%                     |
| Operating           | 5°C ~ 35°C                     |
| Temperature         |                                |
| Operating Humidity  | 30% ~ 85%, R.H. non-condensing |
| Operating Pressure  | 80kPa ~ 106kPa                 |
| Operating Attitude  | 0 – 2,000m                     |

#### Physical Characteristics

| Weight       | 1000g                     |
|--------------|---------------------------|
| Dimensions   | 85(H) x 84(W) x 245(D) mm |
| Package unit | 1ea                       |

#### Accessories

| Basic    | SpO2 Sensor(ACCY-0A0PRB Reusable type) | 1EA |
|----------|--|-----|
|          | Power Cord                             | 1EA |
|          | User Manual                            | 1EA |
| Optional | SpO2 Sensor(pediatric, neonate)        |     |
|          | SpO2 Sensor(pediatric, neonate Y-type) |     |
|          | SpO2 extension cable(9pin)             |     |
|          | IV Pole Clamp                          |     |
|          | Carrying Bag                           |     |

#### Sensor and Cable Length

| Туре                  | Patient Size |
|-----------------------|--------------|
| SpO2 adult Sensor     | 1m, 3m       |
| SpO2 pediatric Sensor | 1m, 3m       |
| SpO2 neonate Sensor   | 1.5m, 3m     |
| SpO2 neonate Sensor   | 1m, 3m       |
| SpO2 extension cable  | 2.5m         |
| Power Cord            | 1.8m         |

#### General

| Display |  |
|---------|--|
| Display | Color TFT Dual LCD                               |
| Size    | 4.3inch(480 x 272 dots), 2.4inch(240 x 320 dots) |

Rev. 2 -51- OP-EN-02

| Brightness             | 1~5   |  |  |
|------------------------|---|--|--|
| LED Lamp               |   |  |  |
| Display power status   | AC Power Status Lamp                          |  |  |
|                        | Battery Power Status Lamp                     |  |  |
| ALARM                  |   |  |  |
| Alarm indicators       | Alarm message,                                |  |  |
|                        | Alarm sound,                                  |  |  |
|                        | Alarm lamp                                    |  |  |
| Alarm level            | High Priority,                                |  |  |
|                        | Medium Priority,                              |  |  |
|                        | Low Priority                                  |  |  |
| Alarm volume           | 0~7   |  |  |
| Alarm duration         | 30, 60, 90, 120 seconds                       |  |  |
| Trends                 |   |  |  |
| Memory                 | Save continuously for 15 days (for 10 seconds |  |  |
|                        | saving period)                                |  |  |
| Display                | Tabular, Graphic                              |  |  |
| External Communication |   |  |  |
| Serial RS-232          | Software Upgrade                              |  |  |

#### **Manufacturer's Declaration**

#### Electromagnetic Compatibility (EMC)

### (Laution

This monitoring system is intended for use by healthcare professionals only. This monitoring system may cause radio interference or may disrupt the operation of nearby equipment, regardless of whether it is CISPR compliant or not. It may be necessary to take mitigation measures, such as reorienting or relocating the monitoring system or shielding the location.

### **Caution**

The use of accessories, Sensors, and cables other than those specified may result in inaccurate reading s of the monitoring system and increased emission and/ or decreased electromagnetic immunity of the monitoring system.

Rev. 2 -52- OP-EN-02

### • Electromagnetic Emissions

| <b>Emissions Test</b>  | Compliance         | Electromagnetic Environment Guidance   |  |  |
|--|--------------------|--|--|--|
| RF emissions<br>CISPR 11                                       | Group 1<br>Class A | Not intended for use in a residential environment. If used in a domestic environment, may not offer adequate protection to radio-frequency communication services. The user may be required to take mitigation measures, such as relocating or re-orienting the equipment. |  |  |
| Harmonic<br>emissions<br>IEC 61000-3-2                         | Class A            |  |  |  |
| Voltage<br>fluctuations/<br>flicker emissions<br>IEC 61000-3-3 | Complies           |  |  |  |

#### Electromagnetic Immunity

| l   | EN/IEC 60601-1-2  | Compliance level  | Electromagnetic   |
|---|---|---|---|
| Immunity test                                       | Test level  | Compliance level  | Environment-Guidance  |
|   |   |   | Floors should be wood,  |
| Electrostatic                                       |   |   | concrete or ceramic tile.   |
| discharge   | ±6 kV contact   | ±6 kV contact   | If floors are covered with  |
| (ESD)   | ±8 kV air   | ±8 kV air   | synthetic material, the relative  |
| IEC 61000-4-2                                       |   |   | humidity should be at least   |
|   |   |   | 30 %.   |
| Electrical fast<br>transient/burst<br>IEC 61000-4-4 | ±2 kV for power supply lines ±1 kV for input/output lines | ±2 kV for power supply lines ±1 kV for input/output lines | Mains power quality should be that of a typical commercial or hospital environment. |
| Surge<br>IEC 61000-4-5                              | ±1 kV line(s) to line(s) ±2 kV line(s) to earth           | ±1 kV line(s) to line(s)<br>±2 kV line(s) to earth        | Mains power quality should be that of a typical commercial or hospital environment. |
| Voltage dips, short                                 | <5 % <i>U</i> T<br>(>95 % dip in <i>U</i> T)              | <5 % <i>U</i> T<br>(>95 % dip in <i>U</i> T)              | Mains power quality should be that of a typical commercial or                       |
| interruptions                                       | for 0,5 cycle   | for 0,5 cycle   | hospital environment. If the  |
| and   | 40 % <i>U</i> T   | 40 % <i>U</i> T   | user of the [ME EQUIPMENT or  |
| voltage   | (60 % dip in <i>U</i> T)                                  | (60 % dip in <i>U</i> T)                                  | ME SYSTEM] requires   |
| variations  | for 5 cycles  | for 5 cycles  | continued operation during  |

Rev. 2 -53- OP-EN-02

| on power       | 70 % <i>U</i> T           | 70 % <i>U</i> T           | power mains interruptions, it is |
|----------------|---------------------------|---------------------------|----------------------------------|
| on power       | 70 70 01                  | 10 70 01                  | power mains interruptions, it is |
| supply         | (30 % dip in <i>U</i> T)  | (30 % dip in <i>U</i> T)  | recommended that the [ME         |
| input lines    | for 25 cycles             | for 25 cycles             | EQUIPMENT or ME SYSTEM]          |
| IEC 61000-4-   | <5 % <i>U</i> T           | <5 % <i>U</i> T           | be powered from an               |
| 11             | (>95 % dip in <i>U</i> T) | (>95 % dip in <i>U</i> T) | uninterruptible power supply or  |
|                | for 5 s                   | for 5 s                   | a battery.                       |
| Power          |                           |                           | Power frequency magnetic         |
| frequency      |                           |                           | fields should be at levels       |
| (50/60 Hz)     | 3 .0A/m                   | 3 .0A/m                   | characteristic of a typical      |
| magnetic field |                           |                           | location in a typical commercial |
| IEC 61000-4-8  |                           |                           | or hospital environment.         |

#### Electromagnetic Immunity Compliance

|               | IEC 60601-1-2     | Compliance level  | Electromagnetic         |
|---------------|-------------------|-------------------|-------------------------|
| Immunity test | Test level        | Compliance level  | Environment-Guidance    |
|               | Frequency of      |                   | Equation for Separation |
|               | Transmitter       |                   | Distance (d)            |
| Conducted RF  | 3 Vrms            | 3 Vrms            | d = 1.2√P               |
| IEC 61000-4-6 | 150 kHz to 80 MHz | 150 kHz to 80 MHz | 150 kHz to 80 MHz       |
|               | 3 V/m             | 3 V/m             | d = 1.2√P               |
|               | 80 MHz to 800     | 80 MHz to 800 MHz | 80 MHz to 800 MHz       |
| Radiated RF   | MHz               |                   |                         |
| IEC 61000-4-3 | 3 V/m             | 3 V/m             | d = 2.3√P               |
|               | 800 MHz to 2,5    | 800 MHz to 2,5    | 800 MHz to 2,5 GHz      |
|               | GHz               | GHz               |                         |

#### • Recommended Separation Distance

| Rated maximum                            | Separation Distance in Meters  |                                |                                 |  |
|--|--------------------------------|--------------------------------|---------------------------------|--|
| output power (P) of transmitter in watts | d = 1.2√P<br>150 kHz to 80 MHz | d = 1.2√P<br>80 MHz to 800 MHz | d = 2.3√P<br>800 MHz to 2,5 GHz |  |
| 0.01                                     | 0.12                           | 0.12                           | 0.23                            |  |
| 0.10                                     | 0.38                           | 0.38                           | 0.73                            |  |
| 1.00                                     | 1.20                           | 1.20                           | 2.30                            |  |
| 10.00                                    | 3.80                           | 3.80                           | 7.30                            |  |
| 100.00                                   | 12.00                          | 12.00                          | 23.00                           |  |

Rev. 2 -54- OP-EN-02

For transmitters rated at a maximum output power not listed above, estimate the separation distance (d) using the equation in the corresponding column, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.



These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.



Portable and mobile RF communications equipment can affect medical electrical equipment.

Such RF equipment should be used no closer to any part of the monitoring system, including cables, than the recommended separation distance calculated from the equation appropriate for the frequency of the transmitter.

#### **Clinical Studies**

#### Overview

This contains data from clinical studies conducted for Neonate Sensor Adult Sensor, Disposable Sensor and Pediatric Sensor used with the ACCURO oximeter.

#### Methods


The current study for Charmcare, Inc. performed 3/18/10 included 10subjects- 7 women and 3 men. The ACCURO oximeter was studied with three different sensors: adult, disposable and neonate. A radial arterial cannula was placed in either the left or right wrist of each subject. Blood gas analysis to determine oxyhemoglobin saturation was performed on an OSM 3® multi- wavelength oximeter (Hemoximeter, Radiometer, Copenhagen, serial 89R0243 N010). This instrument underwent a full factory technician calibration and reference standards check on June 8, 2009. The instrument was also verified to read within specifications with the use of Radiometer Qualicheck standards within 48 hours of the study dates. No subject was anemic (Hemoglobin ≤ 10 gm •dl-1) and only healthy non-smoking individuals of age 21-49 were included in the study.


The study was reviewed and approved by the University of California at San Francisco Committee on Human Research. The approval date was 3/21/09, expiration April 21, 2010. The Approval number is H6301-01706-24. The approval letter is on file at UCSF. Each subject had control data taken at the beginning of each experiment, with two control blood samples drawn while breathing room air. Hypoxia was induced to different levels of oxyhemoglobin saturation (between 70-100%) by having subjects breathe mixtures of nitrogen, room air, and carbon dioxide. Each plateau level


Rev. 2 -55- OP-EN-02


of oxyhemoglobin saturation was maintained for at least 30 seconds and until pulse oximeters readings had stabilized. Two arterial blood samples were then obtained, approximately 30 seconds apart. Each stable plateau therefore was maintained for at least 60 seconds with SpO2 fluctuating by less than 3%. The plateaus were nominally at 100%, room air saturation, 93%, 90%, 87%, 85%, 82%, 80%, 77%, 75% and 70%. A total of 22-24 samples were obtained that the plateaus across this span. Data were recorded by computer. At least 200 data points were collected for each type of oximeter and probe combination studied.


#### Study Results











#### Conclusion

The pooled results indicate that for a saturation range of 60-80% for SpO2, the acceptance criterion was met for the monitoring system when tested with Neonate Sensor Adult Sensor, Disposable Sensor and Pediatric Sensor. The pooled results indicate that for a saturation range of 70-100% for SpO2, the acceptance criterion was met.

Rev. 2 -56- OP-EN-02

## Attachment 1

## **Factory Set**

The default value of the product is as shown below.

#### **Factory Set**

- \* Timing to save
- When altering parameter values: When a user changes each parameter value at her/his discretion it will be saved in memory that after a quick restarting the last saved values will be applied.
- USER CONFIG: Even if the user changes parameter values at his/her discretion validity of this
  change will be preserved for only 30 seconds after powering off the device before restarting it.
  Authorized maintenance personnel can only change USER CONFIG setting(MAINT>USER
  CONFIG>SAVE) that will be always retrieved if restarted more than 30 seconds after any other
  normal users manipulations.

Below shows the parameters that users can changes.

| Menu          | Item              | Factory set     | * Timing to save |
|---------------|-------------------|-----------------|------------------|
|               | PULSE VOLUME      | 4               | When changing    |
|               | FOESE VOLONIE     |                 | parameter value. |
|               | BUTTON VOLUME     | 4               | When changing    |
|               |                   |                 | parameter value. |
|               | BRIGHTNESS        | 3               | When changing    |
|               |                   | 3               | parameter value. |
|               | SCREEN TYPE       | 3               | When changing    |
| CONFIGURATION |                   |                 | parameter value. |
|               | SPO2 AVERAGE(s)   | 8 (neonate :12) | USER CONFIG      |
|               | WAVE FILL         | OFF             | When changing    |
|               |                   |                 | parameter value. |
|               | PLETH BAR         | ON              | When changing    |
|               |                   |                 | parameter value. |
|               | SWEEP SPEED(mm/s) | 12.5            | When changing    |
|               |                   |                 | parameter value. |
|               | Alarm volume      | 4               | USER CONFIG      |
| ALARM         | Alarm duration    | 120             | USER CONFIG      |
|               | SpO2 Alarm        | ON              | USER CONFIG      |

Rev. 2 -57- OP-EN-02

|          | PR Alarm         | ON             | USER CONFIG      |
|----------|------------------|----------------|------------------|
|          | SpO2 High limit  | OFF            | USER CONFIG      |
|          | SpO2 Low Limit   | 90(neonate:85) | USER CONFIG      |
|          | PR High Limit    | Adult : 120    |                  |
|          |                  | Pediatric :160 | USER CONFIG      |
|          |                  | Neonate : 200  |                  |
|          | PR Low Limit     | Adult : 50     |                  |
|          |                  | Pediatric : 75 | USER CONFIG      |
|          |                  | Neonate : 100  |                  |
| PATIENT  | PATIENT          | ADULT          | When changing    |
| FAIILINI | TATIENT          |                | parameter value. |
| TREND    | SAVE INTERVAL(s) | 10             | When changing    |
|          |                  |                | parameter value. |
|          | SHORT TREND(min) | 30             | When changing    |
|          |                  |                | parameter value. |

#### **Maintenance Set**

Below shows the parameters that authorized personnel with password can change.

| Category    | Name           | Factory Default | * Timing to save        |
|-------------|----------------|-----------------|-------------------------|
| MAINTENANCE | ALARM 0 USE    | OFF             | When changing parameter |
|             |                |                 | value. (Authorized      |
|             |                |                 | maintenance personnel)  |
|             | ALARM 0 REMIND | OFF             | When changing parameter |
|             |                |                 | value. (Authorized      |
|             |                |                 | maintenance personnel)  |
|             | LANGUAGE       | ENGLISH         | When changing parameter |
|             |                |                 | value. (Authorized      |
|             |                |                 | maintenance personnel)  |

### **Factory Set By Patient Type**

\* When transferred to Patient mode, the alarm limit value and SPO2 Average value automatically changes to the USER CONFIG set.

| Item             |            | Factory Default |
|------------------|------------|-----------------|
|                  | ADULT      | : OFF           |
| SPO2 HIGH (%)    | PEDIATRIC  | : OFF           |
|                  | NEONATE    | : OFF           |
|                  | ADULT      | : 90            |
| SPO2 LOW (%)     | PEDIATRIC: | : 90            |
|                  | NEONATE    | : 85            |
|                  | ADULT      | : 120           |
| PR HIGH (bpm)    | PEDIATRIC: | : 160           |
|                  | NEONATE    | : 200           |
|                  | ADULT      | : 50            |
| PR LOW (bpm)     | PEDIATRIC  | : 75            |
|                  | NEONATE    | :100            |
|                  | ADULT      | : 8             |
| SpO2 Average (s) | PEDIATRIC  | : 8             |
|                  | NONATE     | :12             |

Rev. 2 -59- OP-EN-02

# Attachment 2 Alarm Message

The level and cause of a single alarm is presented. '\*' mark can be changed according to the value, and when the value does not need to be changed it is set to a fixed value.

#### **Patient Alarm Message**

| Alarm Message | Level  | Cause  |
|---------------|--------|--|
| HIGH SPO2     | Medium | Occurs when the measured value is greater than the high      |
| LOW SPO2      | Medium | limit or greater than the low limit. Check the status of the |
| HIGH PR       | Medium | patient. Also, check if the limit set value is valid.        |
| LOW PR        | Medium |  |

#### **Technical Alarm Message**

| Alarm Message | Level  | Cause  |
|---------------|--------|--|
| SPO2 FAULT    | High   | Not connected with the oxygen saturation measuring part. |
|               |        | The sensor is not connected to the product. Occurs when  |
| SENSOR OFF    | Low    | the sensor is not correctly connected or a problem is    |
|               |        | detected in the sensor.                                  |
| FINGER OFF    | Low    | The sensor is not attached to the patient.               |
| CHECK PROBE   | Low    | Checking the sensor status is required. The light of the |
|               |        | sensor may be exposed to an external light.              |
| LOW BATTERY M | Medium | Battery is low. Connect AC power and charging is         |
|               |        | necessary.   |
| BATT. ERROR   | Low    | Coin battery is low.                                     |
| CLOCK ERROR   | Low    | An error for time value occurs. reset is necessary.      |





# Charmcare Co., Ltd.

(Gasan-dong, Woolim Lions2-cha)714, 2, Gasandigital1-ro, Geumcheon-gu, Seoul, Korea (ZIP 153-787)

Tel.: +82-2-862-5052, Fax:+82-2-862-5065

Home Page: http://www.charmcare.com

User Manual Rev.3 (2015.03.17)

#### EC REP TECNOMED 2000

C/ Palos de la Frontera, 4, 30 28012 - Madrid Spain

Tel: +34-91-530-0117 Fax: +34-91-539-0624